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
CE 1023

Date of the last revision: 04.05.2021

1. Definition and Product Sizes (provided pairwise)

Lipocollect is a moist heat sterilized product, consisting of two glass columns filled with LDL adsorbent. One treatment is performed with 2 products (column pair).

Lipocollect 200, cont. 200 +/- 15 ml of adsorbent, extracorporeal vol. about 180 ml/ column,
Lipocollect 300, cont. 300 +/- 15 ml of adsorbent, extracorporeal vol. about 260 ml/ column,

2. Intended use

Lipocollect columns are intended for multiple lipoprotein apheresis from blood plasma on the same patient. Lipoprotein apheresis is limited to the treatment of severe hypercholesterolemia and dyslipidemia, according to the physician's decision. The treatment is performed by trained professionals.

3. Mode of action

The effect is achieved by adsorptive binding of low density lipoproteins (LDL and Lp(a)) from blood plasma in the extracorporeal circulation. The binding can be reversed by rinsing the columns with hypertonic saline solution during the treatment. Two columns are used simultaneously during treatment and are each loaded and eluted several times. Thus the elimination capacity is high. Other lipoproteins (VLDL and HDL) are bound to a small extent. Albumin, fibrinogen and immunoglobulins are not bound. Cleaning, disinfection or re-sterilisation are not used during the lifetime of the products. The lifetime of the products is 3 years from the date of manufacture.

4. Principle of the procedure

The extracorporeal circulation is established according to the requirements of the primary separation device and anticoagulated with citrate and heparin. The blood plasma is separated from the blood cells by centrifugation or filtration, transferred to the ADAorb® device (Medicap clinic, Ulrichstein) and from there continuously fed to the pre-rinsed columns. There, low-density lipoproteins are adsorbed, whereas the LDL-poor plasma is reunited with the blood cells.

The columns are loaded with rinsing solutions and plasma in a programmed procedure using the ADAorb® device, and are usually loaded and regenerated several times during a treatment. Part of the programmed procedure is the final rinsing of the columns with the color-coded preservative solution, PBS-NaN₃ blue, in which they are stored at 2-5°C.

5. Interactions

Conditioned by the negatively charged adsorbent surface, bradykinin, a physiologic vasodilator, is activated in the plasma inside the columns (see contraindications and side effects). In addition, positively charged drugs that are in the plasma can be adsorbed to the adsorbent and their level of action in the blood is thereby reduced.

6. Contraindications

- Patients with contraindications for the treatment in an extracorporeal circulation;
- Patients with increased bradykinin activity. Such patients are e.g.:
- patients after intake of ACE-inhibitors
- patients with inherited elevated bradykinin activity
- patients with asthma bronchiale
- patients with the rare hereditary angioedema (Quincke oedema);
- Patients with haemorrhagic diathesis as in haemophilia, von Willebrand-Jürgens syndrome, thrombocytopenia or other causes;
- Patients with signs of an infection;
- Patients with a positive infectserology, according to criteria of blood donation

Comment:

Due to insufficient information the products are not suitable for the following cases: the treatment of pregnant women and of children below the age of 16. Patients with implants (such as heart valves or hip prosthesis), transplanted organs or with immune deficiency for bacterial infections may have an increased risk of infection through reusable columns.

7. Side effects

- Known side effects of the extracorporeal circulation are common, especially puncture damage, side effects especially due to hypocalcaemia from the citrate anticoagulant, such as perioral tingling, paraesthesia, chills, nausea, vomiting and tremors, in severe cases cardiac arrhythmias and convulsions. Drop in blood pressure due to volume deprivation may lead to somnolence, fainting, anxiety and dyspnoea.
- The following patients are notably predisposed to drop in blood pressure: first time treated patients, low body weight, beta-receptor blocker treatment and increased sensitivity to bradykinin.
- Bradykinin-mediated side effects appear mainly in the beginning of the treatment, when the first column-treated plasma enters the patient's circulation. At this moment vasodilatation may occur, indicated by flush, tachycardia, headache, and drop in blood pressure and in the most unfavourable case in unconsciousness. In these instances, the extracorporeal blood flow must be interrupted immediately (Do not infuse further activated blood!) until the symptoms have disappeared spontaneously. If necessary, a volume substitution, e.g. 0.9% NaCl solution, must be infused. To prevent this so-called anaphylactoid reaction, take care for the following:
- The intake of ACE inhibitors must be excluded safely before each treatment. 2. The extracorporeal blood flow must be routinely slow during the first column loading until

about 10 min. after the start of the plasma reinfusion. 3. Do not supply diluted plasma to the columns.

- Heparin and the reduction of clotting factors favour bleedings up to 36 hours after the treatment. The blood loss by menstruation may be also increased thereby.
- Iron deficiency by blood loss to the extracorporeal device and by sampling is common in chronic apheresis treatment and usually requires ferritin control and iron supplementation.
- Allergic reactions, such as itching or urticaria may occur.

8. Risks of the treatment

- Situations dangerous to life may arise from gross errors, above all in the handling of the technical equipment, such as confounding or forgetting tubing connections or solutions, bypassing air traps, wrong dosage of the anticoagulants, by inadequate treatment of side effects of the apheresis or by disregard of aseptic handling and storage conditions of the columns.
- Extracorporeal adsorbents in general may be harmful by leaching of soluble or particular matter into the patient`s circulation.
- In patients with an undetected pre-existing bleeding disorder, bleeding may occur after the treatment or may be aggravated. An increased blood loss by the menstruation may be favoured too.
- Removal of any other substances besides LDL by the apheresis may be harmful. These can be medications or endogenous substances that are not known in detail.
- Electrolyte imbalances before the apheresis treatment, especially in liver- and renal insufficiency may result to an elevated risk of hypocalcaemia and hypokalemic alkalosis by citrate anticoagulation and infusion of the alkaline column priming solution (bicarbonate buffer).
- Reusable adsorbents bear the risk of contamination and of an infection of the patient, which may be severe.

9. Application method

General remarks

For new customers and on request Medicollect will conduct user training including the demonstration of the column handling.

The Lipocollect columns must only be used with the prescribed liquids and suitable technical equipment. Furthermore, a suitable infrastructure and personnel experienced in performing apheresis treatments must be available.

Aseptic handling and storage at low temperatures with the colour-coded preservative prevents contamination. It is necessary to wear disposable gloves for all tubing- and fluid connections of the extracorporeal system. To be on the safe side, tests for sterility of the column supernatant (PBS-Azide-blue) according to a recognised procedure (e.g. according to Ph. Eur.) must be carried out on both columns at each treatment.

Each column must be labelled with the patient's surname, first name and date of birth or other unique identification before the first application.

In the event of a reduction in performance or at the end of the shelf life, the products must no longer be used. The maximum application frequency is 52 apheresis/year.

Accessories

1 scalpel for one-time removal of the sealing ring in the metallic Luer connection of the two adsorbers.

The following additional accessories are required and are not supplied by Medicollect:

- The Spectra Optia device for plasmapheresis including accessories. Medicollect has no experience with the use of other plasma separators;
- The ADAorb® unit (Medicap clinic, Ulrichstein) including accessories;
- A column holder from Medicap clinic for fixing the Lipocollect columns to the ADAorb® unit and for the transport to the cold storage;
- A 5 µm disposable filter for filtration of the plasma behind the columns (supplied with ADAorb® disposable tubing set);
- Other materials for apheresis, such as suitable cannulas, three-way stopcocks, connectors, etc.;
- A refrigerator or cold room for storing the columns at 2 - 5° C after opening. The temperature must be monitored.

Solutions for the treatment (not supplied by Medicollect but in the Medicap clinic treatment kit). The solutions are medical devices.

- Priming + Rinsing: Isotonic saline in 5 l bag +1.4% sodium bicarbonate solution in a 3 l bag;
- LDL elution: 5% NaCl solution in 5 l bag;
- Preservation: PBS-0.01% NaN₃, blue solution in a 1 l bottle;
- Anticoagulation: 2 x 250 ml sodium citrate 4% for apheresis; ACDA is less recommended (low pH, glucose content), but can be used as an alternative.

Storage and special steps before the first treatment

Please check visually whether the glass columns are undamaged and the tubes are closed. Columns that arrive open or damaged are non-sterile and must not be used.

Aseptic opening of the columns: In the connection of the metallic Luer connectors there is a sealing ring, which is required during sterilization up to the opening to ensure the sterility of the contents (sterile barrier). To remove the sealing ring, place the adsorber on a clean surface, spray the sealing area with an alcohol-based disinfectant and scratch the ring with the disposable scalpel provided. You can then pull it out easily with a clamp, leaving the Luer connection closed. Do not touch the tubes with the scalpel. Beware of injury, properly dispose the scalpel after use, e.g. in a needle dropping container. After removing the sealing ring, disconnect the upper tube with a stable clamp, unscrew the Luer connection, and attach a sterile stopcock at the inlet and outlet. The stopcocks are used to adjust the liquid level above the adsorbent and for initial rinsing. It must be ensured that not both connections are

open at the same time, otherwise the liquid level above the adsorbent will drop and have to be readjusted.

A separate initial rinse procedure must then be performed once in the life of the columns. It is independent of regular rinsing steps of the first and following treatments.

To do this, attach the column holder to the ADAorb® unit and place the two columns into the holder. Then connect a 5000 ml- 5% NaCl bag with the column inlet by an infusion set and adjust the liquid levels according to Fig. 1. Then flush the contents of the bag halfway through each column, followed by 500 ml of PBS-Azide blue (preservative solution) per column. Alternately you can use the special first rinse program of the ADAorb® device for this. Always seal the columns before storage by welding the disposable tubes with the tube welder provided on the Optia device. Store the columns aseptically sealed at 2-5°C with PBS-NaN₃ blue solution.

Steps for each treatment

During transport and storage, the columns remain in the holder in which they are removed from the cold store before treatment and with which they are attached to the ADAorb® device. The preparation of the devices of the treatment system (Spectra Optia and ADAorb®) is specified by the device manufacturers. The connection of the devices to each other, the column connection to the ADAorb® device and the connection of the rinsing solutions are determined and trained by Medicap clinic. All rinsing steps before, during and after the treatments are defined in the sequence, rinsing volume and rinsing speed in the treatment program (Lipocollect 200 or Lipocollect 300) of the ADAorb® device and can only be changed by the user within defined limits. For treatments, the ADAorb tubing system is connected directly to the metal Luer connectors of the adsorbers without stopcocks, first the upper connection, then the lower connection. It is important to ensure that both connections are not open at the same time, otherwise the liquid level above the adsorbent drops and has to be readjusted. The separation from the tubing system after the treatments is done by welding with a tube welding device.

Anticoagulation

In contrast to a plasma exchange procedure, where most of the citrate is wasted with the removed plasma, in a column procedure all citrate is infused to the patient. This requires a lower citrate ratio (mean 1:30) but more heparin instead (mixed heparin-citrate anticoagulation). If necessary, the proper anticoagulation also must be adapted to the individual patient and controlled by the responsible physician. Therefore, the following recommendation may be unsuitable in individual cases.

The continuous anticoagulation is added to the drawn blood by the plasmaseparator in the preselected ratio. The columns need citrate for a proper function. Heparin alone is not suitable.

Before the start of the extracorporeal circulation, the patient gets an intravenous Heparin bolus of 30-35 IE/ kg body weight, corresponding to about 2500 IE Heparin for a patient of 75 kg body weight. Additionally, continuous heparin application is usually necessary during the treatment. This may be administered by an infusomat or mixed with the citrate anticoagulant. The amount of Heparin mixed to the citrate anticoagulant should result in a concentration of 0,33 IE Heparin/ ml of drawn blood.

Calculation examples for the heparin amount to be added to the citrate anticoagulant:

Example 1

Selected anticoagulant to blood ratio = 1:25 (1+24), volume of the anticoagulant bag 500 ml:
 $0,33 \times 25 \times 500 = 4.125$ IE Heparin have to be added to the bag.

Example 2

Selected anticoagulant to blood ratio = 1:30 (1+29), anticoagulant bag size 250 ml:
 $0,33 \times 30 \times 250 = 2.475$ IE Heparin have to be added to the bag. In this instance, a second anticoagulant bag may contain less Heparin.

In most instances, it is necessary to antagonize the citrate by infusion of calcium (calcium gluconate or calcium chloride). Usually 2,3-4,6 mmol calcium are applied continuously during the treatment by means of an Infusomat and via 3-way stopcock into the return line at the needle. Calcium must never be added to the extracorporeal circulation, because citrate will be inactivated thereby and clotting occurs in the tubing system.

Duration of the treatment:

The mean treatment duration is 2-4 hours.

10. Warnings and precautions

Danger by particle infusion!

- Ensure the correct connection of the columns: The column inlet is at the top and must be Luer female, the column outlet is at the bottom and must be Luer male!
- Never shake or turn round the columns after the column priming or during the treatment!
- Always ensure that the post column plasma is passed through the supplied disposable 5µm-particle filter!
- The columns must not be frozen!
- The columns must not be shaken vigorously!

Danger by elution solution and preserving solution

- The infusion of these solutions may be fatal!
- During the pre-treatment and post-treatment rinsing procedures of the columns, the patient must not be connected to the extracorporeal circuit!

Danger by air embolism!

- Always ensure that the post column plasma is passed through the return line bubble trap of the plasma separator!

Danger by the generation of Bradykinin!

- Exclude intake of ACE-inhibitors before each treatment!
- Use the prescribed solutions in the prescribed sequence and quantity exclusively!
- During the first column cycle, the plasma flow must be slow!
- Do not carry diluted plasma to the columns!

Danger by contamination!

- Do not use the columns if a contamination is suspected according to the appraisal by the physician! Contact the manufacturer!
- Never use damaged or leaking columns! Contact the manufacturer!
- Handle all connections of the extracorporeal circulation with disposable gloves in an aseptic manner!
- After opening, do not store the columns above 5°C!
- Regular checks of the column supernatant for sterility must be performed as described in section 9!

11. Symbols used on the product label

This symbol means „Attention“. It refers to the requirement for the user, to attend the Informations for Use for important safety-related information, which is not indicated on the medical device, such as warnings and precautions.

This symbol means „Attend the Informations for Use!“ . It refers to the necessity for the user, to refer to the Informations for Use.

This symbol indicates the Lot. Each column vessel is a separate Lot.

This symbol means „sterile fluid path“. It indicates that the inward only of products is sterile.

This symbol indicates „Date of production“. It indicates the date, on which the medical device has been produced.

This symbol means „to be used until“.

2020-01-31

This symbol means „brittle, handle with care“. It refers to a medical device, which may break or may be damaged, if handled without care.

This symbol means „manufacturer“. It indicates the manufacturer of the medical device according to the EU guideline (MDR).

Disclaimer: Medicollect supplies products manufactured according to validated manufacturing and sterilisation procedures and intended for reuse according to the product label and instructions for use. Performance and safety have been proven in more than 10 000 treatments.

However, the choice from the available Lipoprotein-apheresis treatment systems on the market is the responsibility of the user and must be based on the individual needs of the patient.

Fig. 1 Adjustment of the fluid level above the adsorbent
Throughout this procedure the bottom stopcock must be closed!